AngelSounds®

DETECTOR DE LATIDOS FETALES

FETAL HEART DETECTOR

Manual de instrucciones

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Name: Fetal Heart Detector

Model: JPD-100S

Manufacturer: Jumper Medical Equipment CO., LTD

Copyright

Copyright 2009. All rights reserved.

Statement:

Jumper Medical Equipment CO., LTD. owns the copyright of this nonpublic instruction book.

All contents in this manual refer to the related product current situation.

Jumper Medical Equipment CO., LTD. reserves the right to adapt the contents of this manual if necessary, without notice.

Jumper Medical Equipment CO., LTD. Reserves the right to the final interpretation of this manual.

"JUMPER" and "ANGELSOUNDS" are trademarks of Jumper Medical Equipment CO., LTD.

Contents

Section 1: Introduction	16
1.1 Overview	16
1.2 Description	16
1.3 Operating principles	16
Section 2: Safety	16
2.1 Indications for use	16
2.2 Contraindications	17
2.3 Domestic use	17
2.4 Safety terms and conditions	17
2.5 Safety alert description	17
2.6 Symbol description	18
Section 3: Operation	19
3.1 Unpacking and inspecting	19
3.2 Setting up the product	19
3.3 Build in battery	20
3.4 Operating knob and indicator light	20
3.5 Preparation	20
3.6 Operation	20
Section 4: Maintenance and after sales service	20
4.1 Maintenance	20
4.2 Recommended maintenance and care	21
4.3 Equipment and accessories cleaning	21
4.4 Cleaning instructions	21
4.5 Disinfection	22
4.6 Recycling the batteries	22
4.7 Authorized repair service	22

Section 5: Specifications and safety	22
5.1 Specifications	22
5.2 Operating mode	23
5.3 Physical dimensions	23
5.4 Environmental requirements	23
Section 6: Accessories	23
6.1 General description	23
6.2 Product accessories	24
Appendix A EMC: Information guidance and manufacturer's declaration.	

Section 1: Introduction

1.1 Overview

Please become familiar with the equipment functions before use.

CAUTION: do not use this product in life supporting situations.

CAUTION: It cannot replace the professional fetal monitor. When the fetal heart rate is abnormal, when the fetal heart cannot be found, or the fetal movement cannot be, the pregnant woman should immediately go to the hospital to seek the doctor's help

1.2 Description

Lightweight, portable detector designed to provide advanced detection of fetal heart rate (FHR) in a wide range of sounds.

Fetus growth and development can be determined through examination of these indices. Valid for daily application of gynecology and clinical departments.

An internal battery powers this device.

1.3 Operating principles

The fetal heartbeat detector consists of a probe (transmitter and receiver) and a signal-processing unit.

The ultrasonic wave is transmitted from a piezoelectric ceramic component at the front of the probe to the uterus of the pregnant woman. The echo is received by another piezoelectric ceramic component at the front of the probe when the ultrasonic wave reaches the fetal heartbeat. Then, it is converted into voltage. This Doppler signal is detected and demodulated from the received signal. The Doppler frequency is consistent with the rhythm of the fetal systole and diastole. Once the cardiac valves vibrate, a Doppler frequency excursion is formed. Then, an output signal of the vibrating valve is transmitted and sent to the loudspeaker as the rhythmic sound of fetal heartbeat.

Section 2: Safety

2.1 Indications for use

The product is normally applied to fetus above 16 weeks growth, difference in pregnant mater. The normal range of fetal heart rate: 110bpm-160bpm.

This product may be applied to the fetus between 9-12 weeks of growth, depending on the pregnancy characteristics.

The operator will be able to hear the FHR through the headphones.
To record the FHR, the detector can be connected to a recorder, which meets the safety requirements described in IEC60601-1.

2.2 Contraindications

No contraindications have been related to this product use.

2.3 Domestic use

Please, consult your doctor before using this product.

2.4 Safety terms and conditions

Please read the warning signs described below:

DANGER: this sign identifies situations that will cause serious physical injury or death.

WARNING: this sign identifies situations that may cause serious physical injury or death.

CAUTION: this sign identifies situations that may cause minor physical injury, product damage, or property damage.

2.5 Safety alert description

DANGER: fire and explosion hazard

Do not operate the device in the presence of flammable gases to prevent possible explosions or fires.

CAUTION: temperature / humidity / high-pressure

Exposing this product to extreme external conditions may compromise operability of the product.

CAUTION: battery disposal

Please recycle or dispose of the batteries according to the federal, state or local laws that apply in your country.

WARNING: use only approved equipment

Do not use batteries, gel, cables, or optional accessories other than those approved by Jumper Medical Equipment CO., LTD.

CAUTION: possible radio frequency (RF) interference.

RF susceptibility from cellular phones, CB radios, and FM two-way radio may cause interference with the product.

WARNING: adjacent and/or stacked equipment

The device should not be used adjacent/stacked to electronic equipment. If necessary, the product should be observed to verify normal functioning.

CAUTION: other equipment regulation

Equipments connected to the device must comply with the respective IEC regulation (i.e.: IEC950 for data processing equipment or IEC 601-1 for medical equipment). In addition, all configurations shall comply with the system standard IEC 601-1-1 regulation. Any person who connects additional equipment to the input/output data port is configuring medical equipment, and therefore, is responsible to ensure that the system meets the IEC 601-1-1 regulation. The product service port is only intended for use during maintenance by authorized service personnel.

CAUTION: exterior cleaning

When disinfecting the device, use a non-oxidizing disinfectant, such as cleaners based on ammonium salts or glutaraldehyde.

CAUTION: operating environment.

The device is designed for indoor use. The operator must ensure that the environment meets the required operating environmental conditions before using.

CAUTION: cold environments

If the device is stored in an environment with a temperature below the operating temperature, the unit should be allowed to warm up to the needed operating temperature before use.

2.6 Symbol description

The following symbols may appear in this manual, on the device, or its accessories. Some of the symbols represent standards and compliances associated with the product and its use.

Consult product and/or accessories instructions.

Warning information

Authorized Representative in the EU

CE Mark: this product meets the Medical Devices 93/42/EEC requirements

Date of manufacture

Manufacturer

Storage Temperature

Humidity

Atmospheric pressure

Upward

Not hooking

Specifies the product serial number

Batch Code

The equipment should be sent to the special agencies according to your local rule for separate collection after its useful life.

Section 3: Operation

This section provides the description for managing the product.

3.1 Unpacking and inspecting

Our intent is that your order delivery is fulfilled in perfect conditions. However, it is recommended to verify the contents of your package against your packing slip.

This product is designed for simple operation and setup. The following items are included in the box:


- 1 Main Unit
- 1 Headset
- 1 Recording Cable
- 1 User Guide

Please check each component when unpacking to ensure that they are in perfect conditions and have not been damaged during transport.

- Check that all components are in the box
- Check that the device and/or accessories have not been damaged and do not have manufacturing defects.
- Contact immediately with customer service in case of damage or manufacturing defects.

3.2 Setting up the product

Controls and indicators

- 1 On/Off and volume knob
- 2 Indicator light
- 3 Headset socket (two)
- 4 Transducer
- 5 Battery compartment cover

3.3 Build in battery

1. Open the battery cover (5) located behind the device.
2. Take the battery connector and plug it to the battery.
3. Close the battery cover introducing first the left side of the cover.

3.4 Operating knob and indicator light

The device incorporates a single knob (1) to turn on, off, and regulate volume. The indicator light (2) shows whether the device is working.

POWER ON

When the device is off, turn the knob (1) clockwise to activate it. The indicator will light up.

POWER OFF

When the device is on, turn the knob (1) to the left to deactivate it. The indicator light will turn off.

VOLUME ADJUSTMENT

To raise the volume, turn the knob (1) clockwise. To decrease the volume, turn the knob (1) to the left.

AUDIO OUTPUT

In the audio sockets, only connect a recorder that complies with the rules of IEC 60601-1

3.5 Preparation

Please follow these guidelines before using the product:

- Turn on the device using the on/off knob.
- Apply coupling gel on the probe faceplate or abdomen.
- Move the transducer slowly over the lower part of the abdomen.

3.6 Operation

Locate the fetus position manually to find out the best direction to the fetal heart. Place the probe faceplate at the best position for detecting the fetal heartbeat ensuring the faceplate is in full contact with the skin. Move the transducer until you get the best audio signal. The fetal heart is usually at 1/3 below the navel line in its earliest stage. As the gestational period increases, it moves upward with a slight deviation to the right or left depending on the fetus. In case of coupling gel absence, you can use water.

Section 4: Maintenance and after sales service

Although maintaining this device is simple, it is recommendable to follow some guidelines to ensure proper operation. This section describes the maintenance and service needed for the device and its accessories.

4.1 Maintenance

WARNING: failure to perform the recommended maintenance guidelines could lead to the product deterioration and possible health hazards. The manufacturer does not, in any manner, assume the responsibility for performing the recommended maintenance guidelines. The sole responsibility rests on individuals, hospitals, or institutions using this product.

The transducer surface is frangible and must be handled with care. Wipe off any remaining gel after use. These precautions will prolong the service life.
Before each use, the operator shall ensure the equipment has no visible damage that could risk patient safety or the equipment operability.

Please, perform the following maintenance to ensure proper functioning:

- Visual inspection
- Proper cleaning of the unit and it's accessories
- Check the battery fuel gauge
- Test the device performance

4.2 Recommended maintenance and care

Carefully inspect the equipment and accessories before installation, once every 12 months and before each service.

- Check the equipment may not have visible physical damage.
- Inspect all external connections for loose or frayed connectors.
- Inspect the graphic displays for marks, scratches or other damages.
- Verify the safety label is clearly legible.

INSTRUCTIONS	INSPECT FOR	RECOMMENDED
Inspect connectors and accessories	Foreign substances	Clean the product and accessories as recommended
	Damage or breakage	Contact our customer service
Examine accessory cables	Foreign substances	Clean the cables as described in Section 5 of this manual
	Broken parts, cracks, damage, or extreme wear, bent connectors and pins	Replace the cable
Examine disposable accessories	Expired Products	Replace by new products

WARNING: if the equipment and/or accessories are damaged, please contact our customer service. If necessary, the product will have to be sent back to us for reparation.

4.3 Equipment and accessories cleaning

Please follow these instructions to clean the product exterior surfaces and batteries:

- Isopropyl alcohol (70% solution in water)
- Mild soap and water
- Sodium hypochlorite (3% solution in water)
- Quaternary ammonium compounds (such as Lysol) (10% solution in water)
- Do not use abrasive cleaners or strong solvents such as acetone-based cleaners.
- Do not mix disinfectant solutions (such as bleach and ammonia).
- Do not clean connectors with bleach.

4.4 Cleaning instructions

- Before cleaning, turn the device off and disconnect from the mains.
- Before applying disinfectant solution, remove all adherent soil (fluids, tissue, etc.), and wipe thoroughly with water.

- Do not immerse the device. Keep the exterior surface clean and free from dirt and dust with a soft cloth soaked in soap solution, and dry immediately with a dry cloth. Remove any remaining gel from the transducer with a soft cloth and neutral soap.
- Soak the cloth before cleaning.
- Do not allow fluids to penetrate the equipment.
- Avoid using cloths or materials that may scratch the equipment.

4.5 Disinfection

Clean the surface of the main unit and the transducer with the above mentioned solutions. Then, wipe the transducer surface with 70% ethanol or alcohol. Dry with a soft cloth.

CAUTION: Do not use low-temperature steam or other high temperature methods to sterilize the device.

4.6 Recycling the batteries

The batteries supported by this product are recyclable. Remove the old battery from the device and follow your local recycling regulations

WARNING: irregular treatment of batteries may cause damage to health and environment.

4.7 Authorized repair service

This product does not incorporate internal components for technical maintenance. Try to solve any technical problem with the troubleshooting table presented above. If you cannot solve the problem, contact the Jumper Medical Equipment Co., LTD customer service department.

NOTE: the warranty will be void upon unauthorized disassembly or service of the equipment.

Section 5: Specifications and safety

This section presents the specifications and safety standards of the product.

5.1 Specifications

NOTE: The following specifications are subject to change and are only a point of reference.

ULTRASOUNDS

Ultrasonic emitting frequency:	3MHZ
Ultrasonic emitting power:	<10mW/cm ²
Overall sensitivity at a distance of 200mm from the transducer plate (Doppler frequency: 300±50Hz, target velocity: 10 cm/s~40cm/s)	≥90dB
Spatial and temporal peak acoustic pressure:	<8.6kPa
Output power	<26mW
Effective area of the ultrasonic transducer active element	5.3cm ² ±20%
Acoustic coupling medium for normal use	ph:5.5~8, Acoustic impedance: 1.5*10 ⁶ ~1.7*10 ⁶ Pa · s/m

AUDIO OUTPUT

Audio output power:	<0.2W
Audio output power:	Ø 3.5mm

BATTERY

Battery voltage:	9V
Type:	IEC6F22 9V alkaline

5.2 Operating mode

Continuous operating

5.3 Physical dimensions

104mm (W) x 120mm (D) x 61mm, Wt: 0.14kg (including batteries)
 W -4.1, D -4.7in, H -2.4 IN, Wt: 31lbs (including batteries)

5.4 Environmental requirements

OPERATING CONDITIONS

Temperature: 5 °C to 40 °C
 Humidity: <80% RH, non-condensing
 Atmospheric pressure: 86kPa to 106kPa

STORAGE AND SHIPPING CONDITIONS

Temperature: -10°C to 55°C
 Humidity: 10% - 95% RH, non-condensing
 Atmospheric pressure: 50kPa to 106kPa

CAUTION: This product is intended for indoor use. The operator must ensure compliance with all environmental requirements described in this manual before using.

CAUTION: If the product is stored in a room with temperature below the operating room, the unit should be allowed to reach the same temperature before using.

Section 6: Accessories

6.1 General description

This section contains a list of accessories and software for the product.

6.2 Product accessories

This product is available in more than twenty languages. For a list of available languages, please contact with your dealer or Jumper Medical Equipment Co., LTD. customer service department.

ACCESSORIES	
Part Number	Description
JP100S-HS13B	Headphones Ø3.5mm
JP100S-RC1.2M	Recording cable Ø3.5mm

Appendix A EMC: Information guidance and manufacturer's declaration.

CAUTION:

The fetal heartbeat detector must follow the guidelines for installation and operation described by the EMC.

EN

Portable and mobile RF communication equipment may affect the functioning of the fetal heartbeat detector.

The fetal heartbeat detector cannot be used together or adjacent to another equipment.

The fetal heartbeat detector is designed for use in electromagnetic environments specified below. The customer or user of the fetal heartbeat detector should ensure that it is used in such an environment.			
Emissions	Test	Compliance	
RF Emissions CISPR11	Group 1	The fetal heartbeat detector uses RF energy only for its internal function. Therefore, its RF emissions are very low and unlikely to cause interference in nearby electronic equipment.	
RF Emissions CISPR11	Class B	The fetal heartbeat detector is suitable for use in all establishments, including domestic establishments, and those directly connected to the public low-voltage power supply network for residential buildings.	
The fetal heartbeat detector is intended for use in the electromagnetic environments specified below. The customer or user of the fetal heartbeat detector should assure that it is used in such an environment.			
Immunity test	IEC test level 60601	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Contact ± 6 kV Aire ± 8 kV	Contact ± 6 kV Aire ± 8 kV	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

The fetal heartbeat detector is intended for use in the electromagnetic environments specified below. The customer or user of the fetal heartbeat detector should ensure that it is used in such an environment.

Immunity test	PIEC test level 60601	Compliance level	Electromagnetic environment guidance
Radiation RF IEC 61000-4-3	3 V/m 80 MHz a 2,5 GHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to the main unit or any part of the fetal heartbeat detector than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P} \text{ 80MHz a 800 MHz}$ $d = 1,2 \sqrt{P} \text{ 800MHz a 2.5 GHz}$ a. where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. b. Interference may occur near the equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guides may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the fetal heartbeat detector is used exceeds the applicable RF compliance level above, the lights should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the fetal heartbeat detector.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable RF communications equipment and mobile and fetal heartbeat detector.

The fetal heartbeat detector is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the fetal heartbeat detector can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the fetal heartbeat detector, as recommended below, according to the maximum output power of the communications equipment.

Maximum rated output power of transmitter W	Separation distance according to the transmitter frequency m		
	150kHz to 80MHz $d=1,2 \sqrt{P}$	80MHz to 800 $d=1,2 \sqrt{P}$	2,5
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	0.73
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Ref.: ASOUNDS-C

Distribuido por/*Distributed by*: Quirumed, S.L.

MedPath GmbH

Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

Shenzhen Jumper Medical Equipment Co., Ltd.

D Building, No. 71, Xintian Road, Fuyong Street,
Baoan, Shenzhen, Guangdong, China

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